



RONGEUR To trim and remove excess alveolar bone after extraction.
To contour alveolar bone after multiple extractions.



RNG-01

ELEVATOR



E11RX E 11RX

E11LX E 11LX

ROOT ELEVATOR



ERE-03 E-03



ERE-01 E-01
Straight



ERE-02 E-02
Curved

EXTRACTION FORCEP Used to remove teeth from alveolar bone

Upper molars



FX10S

10S
Upper Molars



FX53L

53L
1st & 2nd Upper Molars; left
Pointed beak for engaging
bifurcated buccal root.



FX53R

53R
1st & 2nd Upper Molars; right
Pointed beak for engaging
bifurcated buccal root.

Lower molars



FX17

17
1st & 2nd Lower
Molars



FX222

222
3rd Lower Molars

Upper & Lower molars and premolars



FX150

150
Upper Incisors, Canines,
Premolars & Roots,
Universal



FX151

151
Lower Incisors, Canines,
Premolars & Roots;
Universal



Upper anterior

FX1

1
Standard Upper
Incisors & Canines

Upper & Lower premolars of
adult and molars of children



FX101

101
Upper & Lower Premolars, and
Deciduous; Universal

For children



FX150S

150S
Upper Incisors, Premolars
& Roots



FX151S

151S
Lower Primary Teeth
& Roots; Universal

Upper & Lower roots



FX300

300
Upper Roots



FX301

301
Lower Roots


RESTORATIVE

AMALGAM PLUGGER

Pluggers are used to condense filling materials into cavity preparations.
Ends are serrated.



PLG1-2 PLG 1-2



PLG0-1 PLG 0-1

AMALGAM WELL

To hold amalgam before it is placed in
preparation. Made of stainless steel.



DTB-01

DTB-01
Amalgam Well