

LABORATORY

WAXING & CARVING INSTRUMENTS

Durable aluminium handles offer extremely light weight. Color-coded.

L-PKT PKT Set (5pcs)



L-PKT1 PKT 1
Curved & tapered points are used to flow molten wax on.



L-PKT2 PKT 2
Curved & tapered points are used to flow molten wax on.



L-PKT3 PKT 3
Used to perfect and enhance the supplemental and developmental grooves.



L-PKT4 PKT 4
Used to perfect the external contours and remove excess wax at the cavo-surface margins.



L-PKT5 PKT 5
Used to remove excess wax as cusp ridges are developed. Its contour maintains desired convexity at these ridges.



L-LK1 LK 1
Evance Knife No.1



L-LK2 LK 2
Evance Knife No.2



L-LK11 LK 11
Carving knife 11



L-LCV1 LCV 1
Carver 1



L-LCV3 LCV 3
Carver 3



L-LCV7 LCV 7
Carver 7



L-LCV31 LCV 31
Wax Carver



L-LCV32 LCV 32
Roach



L-LCV33 LCV 33



L-LCV34 LCV 34



L-LS1 LS 1
Cement Spatula



L-LS2 LS 2
Porcelane Spatula



L-LS3 LS 3
Wax Spatula

L-RPD RPD
RPD Carver



L-PFM PFM
PFM Carver




L-LCV 3-K LCV 3-K



Brite Sources USA

www.BriteSources.com

1444 Spring Leaf, W. Dundee, IL 60118
Tel: 312-895-1773



Ordering/Pricing: www.BriteSources.com
Email: sales@BriteSources.com

Tel: 312-895-1773
Fax: 847-589-1298