

INSTRUMENT CASSETTE

EFECAL-01

CAL-01

- > with lid
- > electro-polished stainless steel for rust-free & stain-free
- > exterior dimensions : 224x141x25H(mm)
- > 8 instruments capacity (Max. length 222mm)



EFECAN-05

CAN-05

- > with lid
- > electro-polished stainless steel for rust-free & stain-free
- > exterior dimensions : 202x134x22H(mm)
- > 5 instruments capacity & extra space for mallets, retractors etc
- > Can be used as a simple osteotome kit or scaling kit.



EFECAN-12

CAN-12

- > with lid
- > electro-polished stainless steel for rust-free & stain-free
- > exterior dimensions : 268x202x28H(mm)
- > 12 instruments capacity & extra space for mallets, retractors etc



EFECAN-05M

CAN-05M

- > with lid
- > electro-polished stainless steel for rust-free & stain-free
- > exterior dimensions : 190x86x22H(mm)
- > 5 instruments capacity (Max. length 179mm)

DIAGNOSTIC

What is Carpal Tunnel Syndrome?

Carpal tunnel syndrome (CTS) is an irritation of the nerve that enters your wrist and goes to your hand and fingers. Symptoms may include any or all of the following: pain, numbness, tingling, weakness in the thumb second and third fingers; loss of grip strength; shooting pain toward the elbow, shoulder and neck; neck pain and stiffness. The word carpal refers to the 8 bones in the wrist (two rows of four across). The "tunnel" is formed on the top by a ligament which runs across the carpal bones and on the bottom by the carpal bones themselves. A nerve called the median nerve, runs through this tunnel to the first three fingers of the hand. This nerve controls the senses and muscle strength of these fingers.

Not the first, but the best

Since the development of silicone handle for dental hand instruments, it is getting widely used in more and more areas. It seems to have become a new trend now and makes more users thankful to its advantages. Our silicone handle instruments are the best ones. The handles are definitely autoclavable, and they don't slip or wear out in color.

The advantages of silicone handle

1. The cushion on the handle prevents pain and callosity on fingers
2. No-slip surface enables easy-control
3. Effective for preventing carpal tunnel syndrome
4. Color identification is available

